

## **Medina Lab, Penn State University**

*DNA extraction from frozen coral sperm protocol (modified from Sambrook and Russell)*

### **Solutions**

Lysis Buffer (prepare 50 mL):

- 10 mM Tris-Cl (pH=8.0)
- 0.1 M EDTA (pH=8.0)
- 0.5% (w/v) SDS

#### *Lysis Buffer Example (50 mL)*

*Tris-Cl (pH=8) stock 1 M: 500  $\mu$ L*

*EDTA (pH=8) stock 0.5 M: 10 mL*

*SDS stock 10%: 2.5 mL*

*Nuclease free water: to 50 mL*

DNase free RNase (stock: 10 mg/mL)

Proteinase K (stock: 20mg/mL)

Ammonium acetate 7.5 M

EtOH (100%)

Phenol (buffered at pH=8.0)

TE (pH=8.0)

### **Materials**

Frozen coral sperm sample

Mortar and pestle

Liquid nitrogen

## **Protocol**

0. Before start, prepare 10 mL Lysis Buffer with RNase

*Example: to prepare 10 mL Lysis Buffer + RNase*

*add 20  $\mu$ L RNase (stock: 10 mg/mL; final 20 mg/mL) to 10 mL Lysis Buffer*

1. Pulverize coral sperm

- snap-freeze approx. 1 g of coral sperm into stainless-steel mortar filled with liquid N<sub>2</sub>
- grind to fine powder
- allow liquid N<sub>2</sub> to evaporate, then add powder (little by little) to 10 vol of lysis buffer (with RNase) in a 50 mL Falcon tube
- shake the tube to allow powder to submerge
- when all powder is in solution, incubate tube for 1 h at 37°C

2. Gently add proteinase K to final concentration of 100  $\mu$ g/mL

3. Incubate at 50°C overnight (or at least 3-4 h)

4a. Cool solution to room temperature, add equal volume of phenol (pH=8.0).

4b. Gently mix for 10 min (or until solution has formed an emulsion)

5. Spin at 5000 g (6500 rpm in Sorvall SS-34 rotor) at room temperature for 15 min

6. Use wide-bore pipette (3 mm diameter) to transfer viscous aqueous phase to fresh centrifuge tube

7. Re-extract aqueous phase with 1 vol. phenol until solution becomes clear (4-5 times) (Note: may need to precipitate protein and DNA clots by spinning without addition of phenol.)

8. Isolate DNA

- add 0.2666 volumes of 7.5 M ammonium acetate and 2 vol of EtOH (100%), and swirl tube
- remove precipitate using a Shepherd's crook
- wash pellet twice with 70% EtOH
- air-dry and resuspend in 1 mL TE (for 24 to 48 hours)