

## Medina Lab, Penn State University

*Montastraea faveolata* host genotyping  
Last updated 01/22/09 by Shini

\*Status: needs to be updated

### Overview

This protocol describes the genotyping of *M. faveolata* host colonies by using microsatellite markers described in Severance et al. (2004, Molecular Ecology Notes, Vol. 4: 74-76). The general methodology is based on the publication by Schuelke (2000, Nature Biotechnology, Vol. 18: 233-234). In principle, this method can be used to genotype any organism, if microsatellite markers are available. The first step includes the amplification of the target sequence using a primer mix that contains an M13-tagged forward primer, a reverse primer and a fluorescently labeled M13 primer. The PCR products are subsequently analyzed on an ABI3700 sequencer.

### A) PCR-protocol

The method takes advantage of amplifying the target sequence using an M13-tagged forward primer, which is limited in amounts, so that once it is used up, the dye-tagged M13-primer will amplify the M13-tagged target sequence.

*Note: Before you start the amplification of microsatellites with the primer mixes described below, you need to make sure that the microsatellite primers (without the M13-dye tagged primer) produce a single band at the expected size range!*

#### 1. Set up primer mixes

Each primer mix will contain the M13-tagged forward primer (in limited amounts), a reverse primer, and a dye-tagged M13 forward primer.

For each microsatellite, set up:

	Primer stock conc.	1x (2 uL /rxn)*
M13-dye (e.g. FAM,VIC,NED) primer	5 uM	0.6 uL
M13-tagged forward primer	5 uM	0.133333 uL
Reverse primer	5 uM	0.6 uL
Total:		2 uL

\* scale up all primers according to the number of individuals to be genotyped

## 2. Set up PCR reactions

For each reaction, set up:

	1x (20 uL rxn)*
PCR-grade H <sub>2</sub> O	14.4 uL
10 x PCR buffer	2 uL
dNTPs (10 mM)	0.4 uL
Primer Mix (see above)	2 uL
Taq polymerase (5U/uL)	0.2 uL
Template DNA (~20 ng/uL)	1 uL

\* scale up all primers according to the number of individuals to be genotyped and number of microsatellites to be used.

## 3. Run the following program:

Hot-start: 94°C – 2’  
35 cycles: 94°C – 30’  
51°C – 30’’  
72°C – 1’  
1 cycle: 72°C – 2’  
For ever: 4°C

## 4. Check PCR amplicons:

Run out 5 uL of the PCR reaction on a 1% agarose gel.

## 5. Sizing PCR products on ABI 3700

To be updated

## Materials

*M. faveolata* microsatellite primer sequences

maMS8\_M13-fw

CACGACGTTGTAAAACGACTCTTGCCTATCAGCAGAGGAG

maMS8\_rv TCTGCAAACCAATGTACCATCT

maMS11\_M13-fw

CACGACGTTGTAAAACGACCAGACGGATTAAATAGTCTCCCA

maMS11\_rv GACGAATTTTGCCGAGTCAC

maMS12\_M13-fw

CACGACGTTGTAAAACGACGGACCTAAACGGGAACACAA

maMS12\_rv GAAAGGCTATTCAAAGCTGGG

maMS2-4\_M13-fw

CACGACGTTGTAAAACGACTGCTTTGACAGCTACGCAAT

maMS2-4\_rv CCGGGAATTTAGCTATTTGG

maMS2-5\_M13-fw

CACGACGTTGTAAAACGACTTGAAGTAAACAGTACGGAAAGG

maMS2-5\_rv TTCATGTAAACCTGTCGCTGTC

maMS2-8\_M13-fw

CACGACGTTGTAAAACGACCCCCTTTGTCACACATCTTTC

maMS2-8\_rv ATGAAGGATAGGCCGCACT

maMS2-17\_M13-fw

CACGACGTTGTAAAACGACTACCCAGGGCCTCTCTTTTT

maMS2-17\_rv CACGTAATGGCAACGTATGG

M13-Tag:

CACGACGTTGTAAAACGAC