

Medina Lab, Penn State University

RNA Extraction Protocol for Coral Larvae

Last updated: 17 March 2009 by Kevin Portune

Protocol

-Larvae should be stored in RNA later (or some preservative)

1. Label 2 ml screw cap tubes.
2. Fill 2 ml tubes with 0.3 ml each of 0.1 mm glass and 0.5 mm zirconia/silica beads.
3. Pipet off RNA later from larval samples and transfer larvae to 2 ml tubes using a spatula.
4. Add 1 ml Qiazol.
5. Bead-beat for 2 min using Biospec Mini-Beadbeater.
6. Spin at 12,000 x g for 10 minutes at 4°C to pellet larval waste.
7. Transfer supernatant to new tube.
8. Let sit for 5 min at RM temp.
9. Add 1/5 vol (200 µl) chloroform, shake for 15 s and let sit for 3 min at RT°C.
10. Spin for 15 min at 13,000 x g at 4°C.
11. Transfer aqueous phase (~500 µl) to 2 ml RNase-free centrifuge tubes.
12. Add 1 vol. (~500 µl) chloroform, shake for 15 s and let sit for 3 min at RT°C.
13. Spin for 15 min. at 13,000 x g at 4°C.
14. Transfer upper aqueous phase (approx 400 µl) to 2 ml RNase-free centrifuge tubes, add 1 vol. isopropanol (~400-500 µl), and vortex.
15. Incubate for 10 min at RT°C, then spin for 20 min at 13,000 x g at 4°C.
16. Discard supernatant, add 1 ml of EtOH (75%), flick tubes to wash the pellet.
17. Spin for 5 min at 13,000 x g at 4°C.
18. Repeat wash step.
19. Air-dry pellet for a couple of minutes and resuspend in 100 µl RNase-free water.
20. Measure concentration on Nano-drop.
21. Purify total RNA with RNeasy kit (Qiagen).

